

# C-Trak Galaxy

Gamma Probe System




# A perfect combination of sensitivity and directionality for sentinel node biopsies


The C-Trak Galaxy System provides accurate and precise detection of gamma radiation for use during sentinel node biopsies in breast, melanoma and other types of cancer.

## Superior Directionality

The directionality of a gamma probe is not only determined by the collimation, but the energy of the gamma rays. Incoming gamma rays are deflected by scattering within the tissue. This scattered radiation, known as **Compton Scattering**, can give a false representation of where the specimen of interest is, especially if deep within the tissue.



The C-Trak Galaxy system uniquely eliminates scattered radiation, therefore providing superior directionality. The surgeon can then dissect in the right direction to find the tissue of interest, making smaller incisions and dissecting less healthy tissue in the process.



**Energy Threshold and Window Technology**

For a variety of isotopes including Tc-99m and I-125 in radioactive seed localisation.

**Large Touchscreen Display**

Provides maximum visibility.

**Large Count Range**

Up to 100,000 cps for pinpointing high activity sources such as I-125 seeds.

**Fast Calibration and Enhanced Diagnostics**

Quickly ensures the system is functioning correctly.

**Power Outlet Operated.**

**Multiple Probes**

These can be stored within the system for use in a range of procedures.

**Stored Timed Counts**

Stored locally for retrieval of counts at a later date – eliminates the need for manual transcription.

**Built-in User Manual**

Care Wise Medical Products  
C-Trak® Galaxy System

Probe Type: Calibration Probe  
Isotope: Cobalt-57  
Count Time: 10 seconds  
Squelch: Not Set  
Range: 0 - 100

1479  
counts per second

Decrease Range   Increase Range

Take Timed Count   View Timed Counts   Change Count Time   Enable Squelch   Main Menu

C-Trak® Galaxy

# C-Trak Galaxy Gamma Probe System

## Probes and Collimators

### OmniProbe

The OmniProbe has removable collimators making it flexible for a variety of procedures.

The collimator has a 0.2 inch (5 mm) opening which is optimum for the most common procedures such as in breast sentinel node biopsies.

- OmniProbe is available in either straight or angled orientation.
- OmniProbe PET used for detecting F-18 in FDG.
- OmniProbe EL for Laparoscopic use (available in straight, 20° or 90°).

In addition to the OmniProbe the C-Trak Galaxy System can be fitted with alternative probes for specialist applications.

### Lechner Collimator

The additional Lechner collimator has an even smaller field of view 0.1 inch (3 mm) making the probe even more directional but at a compromise of some sensitivity.

The Lechner collimator is especially useful for complex procedures where the nodes are close together and/or close to the injection site such as in the head and neck sentinel node biopsies.



### Specification for Standard OmniProbe

Energy Range	27 - 364 keV
Length / Width	6.7 x 0.6 inches (170 x 15 mm)
Weight	0.14 ounces (140 g)
Sensitivity	Typically 1150 counts/sec/μCi @122 keV (Co-57)
Side Shielding	> 99.9% @140 keV
Probe Tip Diameter	0.5 inch (15 mm) with collimator 0.4 inch (11 mm) without collimator

CsI (Caesium Iodide) detection technology.  
Superior patented collimator technology.

## Optional Accessories

- Printer for efficient record keeping.
- Workstation to store all system components.



### Europe & Worldwide

#### Southern Scientific Limited

Scientific House, The Henfield Business Park  
Shoreham Road, Henfield, BN5 9SL, UK

E-mail: [info@southernscientific.co.uk](mailto:info@southernscientific.co.uk)

Tel: +44 (0)1273 497600

Fax: +44 (0)1273 497626

Web: [www.southernscientific.co.uk](http://www.southernscientific.co.uk)

Version 1.2 February 2018

### USA & Canada

#### C/o LabLogic Systems, Inc.

East Pointe Park, 1040 East Brandon Blvd.  
Brandon, FL 33511-5509, USA

E-mail: [solutions@lablogic.com](mailto:solutions@lablogic.com)

Tel: +1-813-626-6848

Fax: +1-813-620-3708

Web: [www.lablogic.com](http://www.lablogic.com)

A  Group Company