

FLIR GRIFFIN[™] G465

Mobile GC/MS with Real-Time Vapor Survey Mode



The FLIR Griffin G465 is the most flexible transportable GC/MS (Gas Chromatograph/Mass Spectrometer) system available for mobile and forward-deployed labs, reconnaissance vehicles, and other portable platforms. Like the Griffin G460, the G465 features an integrated split/splitless injector and universal sampling port. The addition of a heated sampling line enables survey mode for near real-time vapor threat identification. Industry-leading chemical matching techniques provide high-confidence results in complex environments and eliminate data interpretation in the field. Intuitive user controls, on-screen guidance, and simple data presentation expedite decision making for field operators. The Griffin G465 provides the ease of use and intelligence needed to perform real-time countermeasures that protect public safety.

BEST-IN-CLASS CHEMICAL IDENTIFICATION AND BETTER INFORMED DECISION-MAKING

Guesswork is eliminated with high-confidence chemical matching in complex environments

- Method selection tool guides user through operation
- Automated chemical matching eliminates data interpretation in the field
- Color-coded, go/no-go alarms expedite on-site decision-making

THE MOST FLEXIBLE GC/MS WITH REAL-TIME VAPOR SURVEY MODE

Always mission-ready with unmatched sampling options for vapor, liquid, and solid samples

- Integrated heated sampling line for survey mode missions
- Integrated split/splitless injector accepts traditional syringe, SPME fiber, headspace, and PSI-Probe injections
- Integrated Universal Sampling Port provides continuous air monitoring capability and quick-connect fitting for portable plug-and-play samplers

RUGGED, VEHICLE-MOUNTABLE PLATFORM

GC/MS operates while in motion to expedite mission completion

- MIL-STD 810G design resists shock and vibration
- Internal active vacuum pumping maximizes system uptime
- Integrated shock isolation system minimizes footprint on a vehicle



Vehicle mount with integrated real-time vapor sampling



Universal Sampling Port with quick-connect Griffin X-Sorber[™] portable air sampler

Specifications

Griffin G465	
TECHNOLOGY	
Technology	GC/MS; fully integrated low thermal mass gas chromatograph (LTM-GC) and MS/MS-capable ion trap mass analyzer
Mass Range / Scan Rate	35-425 m/z; up to 10,000 m/z per second @ 20 points per m/z
Ionization Type	Internal electron ionization (EI)
Detector	Conversion dynode electron multiplier
LTM-GC Column	Standard VB-5MS (15 m x .18 mm x .18 µm); others available; programmable 40 to 300°C
Calibrant	Onboard FC-43 (Perfluorotributylamine)
Carrier Gas	Connection for external gas source (choice of He or H ₂); available from many vendors; H ₂ generator available
SAMPLING & ANALYSIS	
Sample Introduction	Split/splitless injector accepts: <ul style="list-style-type: none"> - Direct syringe injection (1 syringe included) - SPME fiber (optional) - Manual headspace sampler (optional) - Autosampler (optional) - PSI-Probe™ thermal separation via TAG™ (optional) - PSI-Probe thermal separation via GERSTEL-Twister® (optional) Universal Sampling Port with standard Tenax TA and Carboxen 1017 dual bed preconcentration tubes accepts: <ul style="list-style-type: none"> - Direct air intake onto precon tubes (adaptor included) - Griffin Purge and Trap for water analysis (optional) - Griffin X-Sorber handheld vapor sampler (optional) Membrane Introduction Mass Spectrometry (MIMS) Inlet: <ul style="list-style-type: none"> - Direct air intake for near real-time air monitoring capability
Sample Phase	Solid, liquid, and vapor
Threats	Detects and identifies explosives, narcotics, CWAs, TICs, environmental pollutants, and other chemicals
Sampling & Analysis	Full identification in 4-15 mins for most chemicals, near real-time in survey mode
SYSTEM INTERFACE	
Display & Alerts	Full automation by connection with computer
Communication	Ethernet connection TCP/IP; remote operation and diagnostics
Data Storage	Data automatically stored on supplied laptop (500 GB)
Simplified User Interface	Griffin System Software (GSS); GriffinLib, NIST and AMDIS mass spectral libraries included; capable of user-defined library
Training Requirements	1-2 days depending on level of training desired; Operator, Developer, and Full System certifications available
POWER	
Input Voltage	100-240 VAC; 24 VDC (+/- 5%, 25 A, 1000 W)
Cold Start Time	<30 mins (includes automatic tuning/calibration)
ENVIRONMENTAL	
Operating Temp / Humidity	41 to 104°F (5 to 40°C); <85% relative humidity
Storage Temp	-13 to 131°F (-25 to 55°C)
PHYSICAL FEATURES	
Dimensions (L x W x H)	19.2 x 19.2 x 21.1 in (48.8 x 48.8 x 53.6 cm)
Weight	99.5 lbs (45.0 kg)
Enclosure & Protection	Rugged, internal shock mounting system; integrated vacuum system contains mini turbomolecular pump & quad diaphragm; no external shock table or vacuum system required



HEADQUARTERS
 FLIR Systems, Inc.
 27700 SW Parkway Ave
 Wilsonville, OR 97070

DETECTION SALES, AMERICAS
 FLIR Detection, Inc.
 2800 Crystal Drive, #330
 Arlington, VA 22202
 Phone: +1-877-692-2120
 detection@flir.com

DETECTION SALES, APAC
 FLIR Detection, Inc.
 3 Pickering Street #03-49
 Nankin Row
 Singapore - 048660
 Phone: +65-6822-1596
 detection@flir.com

DETECTION SALES, EMEA
 FLIR Detection, Inc.
 Luxemburgstraat 2
 2321 Meer
 Belgium
 Phone: +32 (0) 3665 5106
 detection@flir.com

www.flir.com
 NASDAQ: FLIR

Equipment described herein may require US Government authorization for export purposes. Diversion contrary to US law is prohibited. Imagery for illustration purposes only. Specifications are subject to change without notice. ©2016 FLIR Systems, Inc. All rights reserved. 16-1119