



Coriolis® RECON



Portable air sampler for biological warfare agents detection

- High air flow rate for collection of all pathogens > 0.5 μ m
- Liquid sample, compatible with any type of analysis
- Ruggedized equipment for harsh environmental conditions
- Portable, easy-to-deploy and to set-up



WWW.BERTIN-INSTRUMENTS.COM





Coriolis® RECON

A ruggedized bio-aerosol sampler


Coriolis® RECON is a ruggedized bio-aerosol sampler, dedicated to CBRN Recon teams or first responders, with quick deployment in case of bio-threat suspicion. Coriolis[®] RECON is efficient, portable and has been ruggedized for unfamiliar environment. Coriolis[®] RECON offers flexibility with a choice of operating modes:

• Autonomous sampling: triggered by the operator wearing IPE; Coriolis® RECON can be used by Recon teams or first responders or for mobile applications to rapidly obtain a sample and identify the biological threat.

- Biological sentry mode: Coriolis[®] RECON can be set up in a standby mode, waiting for an order of a warning system.
- Long time surveillance: Coriolis[®] RECON can sample air during up to 6 hours (Long Time Monitoring option).


Long Time Monitoring Option is available for the Coriolis[®] RECON. This option is ideal for the monitoring of a special area like Remote Base Surveillance. This option adjusts automatically the volume of collection liquid for long sampling during up to 6 hours. The volume of liquid inside the cone is always optimal for an efficient sampling.

Wet Cyclone - a patented technology



Collection liquid is injected into the cone while aspirating.

Air is drawn in, thanks to a whirling motion to form a vortex.



Particles are pushed against the wall due to centrifugal force, and separated from air to be concentrated into the liquid.



CBRN Products Range

Discover the other products of Bertin Instruments





GammaTRACER Spider Quick deployment probe



APPLICATION

PRINCIPLE Wet cyclone

COLLECTED PARTICLES SIZE > 0.5 µm

AIR FLOW RATE 600 L/min

COLLECTION TIME

WATERTIGHT Yes – IP<u>5</u>4

LIQUID SAMPLE 20 ± 5 ml

DIMENSIONS (L X W X H) 365 x 220 x 306 mm without cane

WEIGHT 10 kg (20 kg with case)

AUTONOMY

MAIN SUPPLY 100 - 240 V

OPERATING TEMPERATURE pack option

Bertin Technologies S.A.S • Parc d'Activités du Pas du Lac • 10 bis, avenue Ampère, 78180 Montigny Le Bretonneux WWW.BERTIN-INSTRUMENTS.COM

RCS: 422 511 204 VERSAILLES | SIREN: 422 511 204 | Lessee manager of Saphymo's on-going business