



EASYBOX

SHIELDED CARRIER
FOR SYRINGES OF
RADIOPHARMACEUTICAL SUBSTANCES



Robustness: shock resistance.

Ergonomics: optimised weight, swivelling handle for easy opening, lockable to prevent swinging during transport.

Easy decontamination and cleaning thanks to smooth surfaces made of strong plastic.

Characteristics

General

Radiation protection:
3 mm lead on 4 sides and
6 mm lead on the side ends

Overall dimensions:
L 255 x D 142 x H 145 mm

Internal dimensions:
L 195 x D 100 x H 40 mm

Weight: 3.3 kg.

Package

Package dimensions:
L 320 x D 240 x H 170 mm

Package weight: 5 kg.

Ref.: .00013222



Download
the product data sheet
lemerpax.com