



PFE / PME VIAL SHIELD

LEAD GLASS VIAL SHIELD



Optimum radiation protection: fully lead glass and tungsten design for minimum operator exposure.

Comfort: content visible at 360°.

Compatibility: adapts to all types of vial (through the use of a plastic adaptor).

Characteristics

General	PFE	PME
Overall dimensions:	Ø 57 x H 115 mm	Ø 67 x H 115 mm
Internal dimensions:	Ø 31 x H 62 mm	
Radiation protection (mm lead):	10	15
Weight:	1.4 kg.	1.9 kg.
Package		
Package dimensions:	L 210 x D 150 x H 150 mm	
Package weight:	1.5 kg.	2 kg.
Ref.:	.02250010	.02250011



Download
the product data sheet
lemerpax.com



PHE VIAL SHIELD

LEAD GLASS VIAL SHIELD



Optimum radiation protection: fully lead glass design for minimum operator exposure.

Comfort: content visible at 360°.

Compatibility: adapts to all types of vial (through the use of a plastic adaptor).

Characteristics

General	PHE
Overall dimensions:	Ø 86 x H 115 mm
Internal dimensions:	Ø 80 x H 62 mm
Radiation protection (mm lead):	25
Weight	2,63 kg
Package	
Package dimensions:	L 320 x P 240 x H 170 mm
Package weight:	3,5 kg
Réf. :	02250012



Download the product data sheet
lemerpax.com