



EASYVIEW

COMPLETE RANGE OF ERGONOMIC TUNGSTEN SYRINGE SHIELDS

Radiation protection :

- Tungsten thickness: min. 2 mm
- Lead glass thickness: 8 mm (density 5,2)
- Mitigation >99%

Range compatible with all syringes on the market.

Ergonomics: easy gripping and extra large viewing window with zoom effect.

Easy maintenance: quick and easy glass changing by removing a single screw.



PREPARATION

INJECTION

STORAGE

PRODUCTS +

Characteristics

General	1 mL			2 mL		2.5 mL	3 mL	5 mL			10 mL			
Length:	70 mm	75 mm	76 mm	50 mm	53 mm	59 mm	69 mm	57 mm	61 mm	64 mm	73 mm	75 mm	83 mm	85 mm
Inner dia.:	9 mm	11.5 mm	9 mm	11.5 mm			15 mm			18 mm				
Radiation protection:	2 mm tungsten / 8 mm lead glass (density 5,2)													
Package dimensions:	L 130 x D 100 x H 50 mm													
Package weight	0.2 kg + weight of syringe shield													
Ref. :	00021896	00022655	00022638	00021314	00022641	00021318	00022652	00022311	00022631	00022644	00021289	00022658	00022661	00027962



Download the product data sheet
lemerpax.com



EASYVIEW HE

COMPLETE RANGE OF ERGONOMIC TUNGSTEN SYRINGE SHIELDS



Compatibility: range compatible with all syringes on the market.

Non-magnetic syringe shield adapted to PET MRI.

Ergonomics: easy gripping.

Easy maintenance: quick and easy glass changing by removing a single screw.

PREPARATION

INJECTION

STORAGE

PRODUCTS +

Characteristics

General	1 mL			2 mL		2.5 mL	3 mL	5 mL			10 mL			
Length:	70 mm	75 mm	76 mm	50 mm	53 mm	59 mm	69 mm	57 mm	61 mm	64 mm	73 mm	75 mm	83 mm	85 mm
Inner dia.:	9 mm	11.5mm	9 mm	11.5 mm				15 mm			18 mm			
Radiation protection:	6 mm tungsten / 8 mm lead glass													
Package dimensions:	L 130 x D 100 x H 50 mm													
Package weight	0.2 kg + weight of syringe shield													
Ref. :	00024946	00024953	00024947	00024948	00024949	00024950	00024952	00024954	00024951	00024955	00024956	00024957	00024958	00028075



Download the product data sheet lemerpax.com