



MODEL 8188



MODEL 8188 EXTENDED


Krypton 81m is most efficiently delivered using a reservoir system (JNM 1983: 24). This NP Krypton 81m Delivery System has been designed by Peter Childs and Neil Franklin, of City Hospital, Birmingham.

In order to remove the need for any tubing support such as a drip stand, the reservoir has been reduced to the minimum.

These easy to use, single use systems ensure quick efficient krypton delivery to the patient. The standard system (pictured above) is supplied in boxes of 50. Each box includes 10 thin bore delivery tubes. Also available are complete systems with a mouthpiece and nose-clip, or with a medium adult mask.



MODEL 8189



MODEL 8191

NP KRYPTON 81M DELIVERY CIRCUIT ACCESSORIES




If the kits do not reflect your own preferences, these accessories can be purchased separately and may be used with the standard circuit, model 8188.

For further information on prices or to arrange a demonstration, please contact Rita Coleman on:

01273 497600

5070	Large Adult Mask	Box of 20
5071	Adult Mask	Box of 20
5072	Small Adult Mask	Box of 20
0142	TRU-FIT Mouthpiece and Nose-clip	Box of 20
1604	Disposable Mask Harness	Box of 20
LC2	Air Pump (as recommended by supplier of Krypton Generators)	Box of 1
	Remote switch for Air Pump	Box of 1



Scientific House, The Henfield Business Park,
Shoreham Road, Henfield, West Sussex, BN5 9SL
Tel: +44 (0)1273 497600
Email: info@southernscientific.co.uk