



EASYBOX HE CARRIER

SHIELDED CARRIER
FOR SYRINGES OF
RADIOPHARMACEUTICAL SUBSTANCES



PREPARATION

INJECTION

STORAGE

PRODUCTS +

Robustness: shockproof.

Easy hygiene: Easy decontamination and cleaning thanks to smooth surfaces made of strong plastic.

Ergonomics: optimised weight, swivelling handle for easy opening, lockable to prevent swinging during transport.

Characteristics

General

Overall dimensions:
L 255 x D 142 x H 145 mm

Internal dimensions:
L 185 x D 96 x H 40 mm

Weight: 5.7 kg

Radiation protection:
6 mm lead on 4 sides and
12 mm lead on the side ends.

100 % recyclable

Package

Package dimensions:
L 320 x D 240 x H 170 mm

Weight: 7.4 kg

Ref.: 00014236



Download
the product data sheet
lemerpax.com