



The
XPERT[®] 80-L
Cabinet X-ray System

Designed for large sample imaging


Applications Include:




Life Sciences



Agriculture



Non-Destructive Testing



Forensic Science

High Resolution X-ray System with our Largest Imaging Area Available.

Digital detector

Sizes available up to 17" x 17" (43 x 43 cm).

Microfocus X-ray imaging

90 kV source for maximum clarity and detail.

Integrated HD optical camera

For side-by-side comparisons with X-ray images.

The Image Blender™

Combines X-ray and optical image for a unique view of your region of interest.

Simple to use, fully shielded

No specialized training required. Simple one button operation. Benchtoper model available.

Specifications:

Size (W x D x H)	32" x 37.5" x 59" or 71" (81 x 95 x 150 or 180 cm)
Optical Camera	Included
Detector Size	up to 17" x 17" (43 x 43 cm)
Image Blender	Included
Magnification	Up to 3.5x

