

# The XCELL® 180 Benchtop Irradiator System



# IT'S TIME FOR A 180

On Irradiation Technology

## **Revolutionize Your Research**

with the benchtop irradiator that not only gives you 180 kV, the highest energy available, but also a comprehensive range of innovative proprietary technologies to help you to obtain high quality scientific results and improve efficiency.

### **Features:**

#### **Unmatched Energy Range**

Up to 180 kV for faster and more effective irradiation.

#### Live Feed

Enables an interior real-time view of the chamber to ensure sample placement is correct. Provides confidence that your work is progressing as planned.

#### **Secure Mobile App**

Peace of mind at your fingertips. Check on the status of your system and research from anywhere.

#### **Multiple User Login**

Allows the system admin to easily manage users, and access master event logs.

#### **Integrated Turntable**

Ensures a uniform dose at all intensity levels during the irradiation process.

#### Sleek Design with Easy Setup

No special training required, fast onboarding process for the system - just plug n' play.

All of these features are packed within a small footprint - 22" (W) X 27" (D) X 36" (H) - that still affords an ample homogeneous irradiation area.

CIRP-eligible in the USA

