



# MEDICLIC

## RANGE OF TUNGSTEN SYRINGE SHIELDS FOR LOW AND MEDIUM ENERGY RADIOPHARMACEUTICAL SUBSTANCES

**Compatibility:** with all syringes on the market.

**Simplicity:** can be completely disassembled by the operator for decontamination or to change the glass.

**Ergonomics:** Lead glass window with zoom effect and white lining for optimum viewing.

**Rapidity:** quick and easy syringe insertion/withdrawal system (1 click).



### Consumables

Mediclic tungsten glass  
TC 1 mL short and long  
Ref. 11452

Mediclic tungsten glass  
TC 3 mL long  
Ref. 11492

Mediclic tungsten glass  
TC 2 mL and 2 mL "BD  
Emerald"  
Ref. 11464

Mediclic tungsten glass  
TC 5 mL short and long  
Ref. 11472

Mediclic tungsten glass  
TC 2.5 mL  
Ref. 11462

Mediclic tungsten glass  
TC 10 mL short  
Ref. 11484

Mediclic tungsten glass  
TC 3 mL short  
Ref. 11494

Mediclic tungsten glass  
TC 10 mL long  
Ref. 11482

### Characteristics

General	1 mL short	1 mL long	2 mL	2 mL "BD Emerald"	2.5 mL	3 mL short	3 mL long	5 mL short	5 mL long	10 mL short	10 mL long
Size planned for the syringe:	70 mm	74 mm	52 mm	49 mm	58 mm	62 mm	68 mm	59 mm	62 mm	73 mm	82 mm
Radiation protection:	2 mm tungsten / 6 mm lead glass										
Package dimensions:	L 240 x D 180 x H 50 mm										
Package weight:	0.130 kg + weight of syringe shield										
Ref.:	11453	11451	11463	11465	11461	11493	11491	11473	11471	11483	11481



Download  
the product data sheet  
[lemerpax.com](http://lemerpax.com)



# MEDICLIC HE

## RANGE OF TUNGSTEN SYRINGE SHIELDS FOR RADIOPHARMACEUTICAL SUBSTANCES

**Compatibility:** with all syringes on the market.

**Simplicity:** can be completely disassembled by the operator for decontamination or to change the glass.

**Ergonomics:** lead glass window with zoom effect and white lining for optimum viewing.

**Rapidity:** quick and easy syringe insertion/withdrawal system (1 click).



### Characteristics

General	1 mL short	1 mL long	2 mL	2 mL "BD Emerald"	2.5 mL	3 mL short	3 mL long	5 mL short	5 mL long	10 mL short	10 mL long
Size planned for the syringe	70 mm	74 mm	52 mm	49 mm	58 mm	62 mm	68 mm	59 mm	62 mm	73 mm	82 mm
Radiation protection:	5 mm tungsten / 9 mm lead glass										
Package dimensions:	L 240 x D 180 x H 50 mm										
Package weight:	0.130 kg + weight of syringe shield										
Ref.:	11653	11651	11663	11666	11661	11693	11691	11673	11671	11683	11681

### Consumables

Mediclic tungsten glass TEP 1 mL short and long

Ref. 11652

Mediclic tungsten glass TEP 2 mL and 2 mL "BD Emerald"

Ref. 11664

Mediclic tungsten glass TEP 2.5 mL

Ref. 11662

Mediclic tungsten glass TEP 3 mL short

Ref. 11694

Mediclic tungsten glass TEP 3 mL long

Ref. 11692

Mediclic tungsten glass TEP 5 mL short and long

Ref. 11672

Mediclic tungsten glass TEP 10 mL short

Ref. 11684

Mediclic tungsten glass TC 10 mL long

Ref. 11682

